

SOLVFIL™

Filter Degasser for HPLC

- Save Time Filtering and Degassing HPLC Mobile Phases
- Improve Quality of HPLC Data
- BPA Free Polypropylene Disposable Bottle Top Filter
- Made for HPLC Solvent Chemical Compatibility



1 Liter Size



500ml Size

For the proper care of HPLC Columns, instruments & reproducible methods, it is a “best practice” to filter & degas HPLC mobile phase before use. SOLVFil™ makes this process easy and convenient by offering a 90mm diameter membrane surface area and large BPA free polypropylene reservoirs for faster & better filtering than with traditional glass filtering devices.

SOLVFil™ very quickly clips on to the SOLVFIL™ Cradle Ring or SOLVFil™ Pedestal or uses the included vacuum tube hose barb for connection to vacuum. This filter degasser also directly screws on to any GL-45 threaded mobile phase bottle eliminating decanting from a standard filtration flask into a reservoir bottle; prevents spilling and re-gassing of dissolved oxygen.

For more product information & specifications visit www.cornerstonescientific.com





1L Funnel



500ml Funnel

SOLVFiL™ bottle top filters are shown in a SOLVFiL™ Cradle-Ring on a standard stand (SOLVFiL™ Pedestal is not shown). SOLVFiL™ is compatible with 100ml to 20L reservoir bottles

Product Code	Description	Qty.
CS10-325-5441-5	SOLVFiL™ 1000ml - Nylon 0.45µm HPLC Solvent Filter	12/CS
CS10-325-1441-1	SOLVFiL™ 500ml - Nylon 0.45µm HPLC Solvent Filter	12/CS
CS10-325-5241-5	SOLVFiL™ 1000ml - Nylon 0.2µm HPLC Solvent Filter	12/CS
CS10-325-1241-1	SOLVFiL™ 500ml - Nylon 0.2µm HPLC Solvent Filter	12/CS
CS10-325-5451-5	SOLVFiL™ 1000ml - PTFE 0.45µm HPLC solvent Filter	12/CS
CS10-325-1451-1	SOLVFiL™ 500ml - PTFE 0.45µm HPLC Solvent Filter	12/CS
CS10-1122-RLS-1	Cradle Ring for Bottle Top Vacuum Filtration	1 E/A
CS10-1116-RLS-1	Pedestal for Bottle Top Vacuum Filtration	1 E/A

SOLVFiL™ is a trademark of Chrom-Supply Co.

