

SafeWASTE™

Systems for HPLC With SafeVENT™

Provide a safer work environment by minimizing organic solvent spills in the lab, during disposal and by trapping vapors (VOC) from HPLC solvents

- Proper Care and Disposal of Hazardous HPLC Waste Materials
- Contain VOC - Volatile Organic Compounds from up to 6 HPLC's
- Heavy Duty or UN/DOT Certified Carboys & Containers with Vent-Filters



UN/DOT SafeWASTE™ Container with SafeVENT™

Certified UN/DOT SafeWASTE™ systems are configured for single use by disposal of the container with the supplied closed cap when full. The system contains one BPA free PE container with one 70mm Swivel Cap Kit which includes cap with turning wings using an independent center tubing holder that has 6 1/4 - 28 ports and fittings to accommodate either 1/8th or 1/16th tubing from your HPLC instrument. Solid plugs for ports not in use (in order to maintain a closed system), one SafeVENT™ charcoal filter with one 90° elbow fitting and one single use 70mm solid cap for disposal are included.

The Heavy Duty (HD) SafeWASTE™ system is designed for pouring collected liquids into a larger central waste container for disposal when full. This system contains one BPA free, easy to grip, ergonomic carboy with one 83mm swivel cap kit which includes cap with turning wings, an independent center tubing holder that has 6 1/4 - 28 ports and fittings to accommodate either 1/8th or 1/16th tubing from your HPLC instrument and a port for the vent, 6 each solid plugs for ports not in use (in order to maintain a closed system), one SafeVENT™ charcoal filter and one 90° elbow fitting along with one 83mm solid cap for disposal



HD SafeWASTE™ carboy with SafeVENT™

PWWroduct Code	Description
CM10-332-3411-1	10 Liter HD SafeWASTE™ System
CM10-332-AB01-0	13.5 Liter UN/DOT SafeWASTE™ System
CM10-332-4411-1	20 Liter SafeWASTE™ System

For more product information & specifications visit www.cornerstonescientific.com

